

DRILL PATH VISUALIZATION

SLEEVE LABELING	ADJUSTMENT DRILL PATH-Ø(MM)	COMMENTS
Astra Tech		
Facilitate™ - small	7.00	For implants with Ø 3.5 mm & 4.0 mm
Facilitate™ - large	7.50	For implants with Ø 4.5 mm & 5.0 mm
BIOMET 3i™		
Navigator® - small	6.20	For implants with max. Ø 3.25 mm & 4.00 mm
Navigator® - large	7.20	For implants with max. Ø 5.00 mm
CAMLOG®		
CAMLOG® Guide System	7.00	For implant Screw Line Ø 3.8 mm & 4.3 mm
Dentaurum		
tioLogic pOosition® - Ø 3.3 mm	6.10	
tioLogic pOosition® - Ø 3.7 mm	6.10	
tioLogic pOosition® - Ø 4.2 mm	6.60	
tioLogic pOosition® - Ø 4.8 mm	7.20	
tioLogic pOosition® - Ø 5.5 mm	7.90	
Friadent		
Expertease Xive - ND	6.50	For implants with Ø 3.0, 3.4 & 3.8 mm
Expertease Xive - WD	7.00	For implants with Ø 4.5 mm
Expertease ANKYLOS®	6.6	
MIS® Implant Technologies		
MGUIDE	7.00	

DRILL PATH VISUALIZATION

SLEEVE LABELING	ADJUSTMENT DRILL PATH-Ø(MM)	COMMENTS
Nobel Biocare™		
NobelGuide™ - NP	5.70	For implants with NP platform
NobelGuide™ - RP	7.00	For implants with RP platform
NobelGuide™ - WP	8.00	For implants with WP platform
SICAT		
SICAT Pilothülsen	4.00	
SICAT Hülse-in-Hülse System	5.50	
Straumann®		
Straumann® Guide - Ø 5.0 mm	7.30	
Straumann® Guide - Ø 2.8 mm	4.80	If Ø 5.0 mm sleeve does not fit
ZimmerDental		
Zimmer® Tube A	6.41	For implants with Ø 3.0 mm & 3.7 mm
Zimmer® Tube B	7.50	For implants with Ø 4.1 mm & 4.7 mm